

Product Datasheet

Lysergic acid diethylamide (LSD) Antibody 100911

Product Name	Lysergic acid diethylamide (LSD) Antibody
Catalog Number	100911
Description	Protein A purified monoclonal mouse antibody, cultured in vitro under conditions free from animal-derived components.
Species/Host	Mouse
Tested Applications	ELISA, LF
Immunogen	LSD-BgG
Brand	Medix Biochemica
Form/Appearance	Liquid, may turn slightly opaque during storage
Concentration	≥ 1.5 mg/ml
Storage	+2-8 °C
Note	This product was previously available with catalog code P01-99-66M-P. It has been produced using the same original mouse monoclonal hybridoma cell line.
Isotype	IgG1
Clonality	Monoclonal
Purity	≥ 95 %
Buffer	0.05 M Potassium Phosphate Buffer, pH 7.2 with 0.85 % NaCl, 0.05 % NaN ₃ as a preservative
Specificity	Antibody recognizes lysergic acid diethylamide (LSD)
Matched Pairs	LSD-BSA Antigen LA783
Shelf Life	24 months