

Product Datasheet

Serum, Direct Bilirubin Samples 0.1 - 0.4 mg/dL mg/dL 991-24-PS-DBIL2

Product Name	Serum, Direct Bilirubin Samples 0.1 - 0.4 mg/dL mg/dL
Catalog Number	991-24-PS-DBIL2
Brand	Lee Biosolutions
Form/Appearance	Liquid
Storage	-20°C
Source	Single Human Donor
Shelf Life	Due to differing stabilities of various biomarkers in this complex biological material it is not possible to assign a single expiration date that is relevant to all biomarkers. The end user should monitor the biomarkers of interest to verify suitability over time.