

# Product Datasheet

## Serum, Single Donor, Direct Bilirubin 0.1-0.4 mg/dL 991-24-PS-DBIL2

Product Name	Serum, Single Donor, Direct Bilirubin 0.1-0.4 mg/dL
Catalog Number	991-24-PS-DBIL2
Description	Each sample is a minimum of 1mL per donor. Inquire for larger volumes per donor or custom criteria.
Tested Applications	QC Materials
Alternative Names	Serum, fraction, clotted, post-coagulation, coagulated, powder, sera, blood serum, defibrinated, direct, conjugated, bilirubin, liver
Brand	Lee Biosolutions
Form/Appearance	Liquid
Concentration	Direct Bilirubin, 0.1-0.4 mg/dL
Storage	-20°C
Source	Single Human Donor
Shelf Life	Due to differing stabilities of various biomarkers in this complex biological material, it is not possible to assign a single expiration date that is relevant to all biomarkers. The end user should monitor the biomarkers of interest to verify suitability over time.
Preservatives	None added, inquire if needed