
Product Datasheet

BioCon II D - Defibrinated, Delipidated Plasma H2700

| | |
|----------------------------|--|
| Product Name | BioCon II D - Defibrinated, Delipidated Plasma |
| Catalog Number | H2700 |
| Species/Host | Human |
| Tested Applications | QC Materials, Base Matrix |
| Brand | Bioresource Technology |
| Form/Appearance | Liquid |
| Storage | -20 °C |
| Purity | 0.2 micron filtered plasma |
| Source | Pooled Human Donors |