

Product Datasheet

Anti-Aspergillus V5145

Product Name	Anti-Aspergillus
Catalog Number	V5145
Description	Purified from culture medium by protein A chromatography
Species/Host	Rabbit
Tested Applications	ELISA, IHC, IFA
Brand	ViroStat
Form/Appearance	Purified from culture medium by protein A chromatography
Concentration	≥ 1 mg/ml
Storage	+2-8 °C or -20 °C. Avoid freeze-thaws. Maximum one freeze-thaw during storage.
Isotype	IgG
Clonality	Monoclonal
Purity	≥ 90 %
Buffer	0.01M Phosphate Buffered Saline, pH 7.5, containing 0.1% sodium azide
Cross Reactivity	Reactive with 9 aspergillus species including 6 strains of A. fumigatus. A more detailed table of cross-reactivity is available upon request
Specificity	Secreted hyphal antigen from all Aspergillus species tested
Matched Pairs	Self pairs in ELISA
Shelf Life	4 years