

## Product Datasheet

### Anti-Aspergillus V5145

<b>Product Name</b>	Anti-Aspergillus
<b>Catalog Number</b>	V5145
<b>Description</b>	Purified from culture medium by protein A chromatography
<b>Species/Host</b>	Rabbit
<b>Tested Applications</b>	IHC, ELISA, IFA
<b>Brand</b>	ViroStat
<b>Form/Appearance</b>	Purified from culture medium by protein A chromatography
<b>Concentration</b>	≥ 1 mg/ml
<b>Storage</b>	+2-8 °C or -20 °C. Avoid freeze-thaws. Maximum one freeze-thaw during storage.
<b>Isotype</b>	IgG
<b>Clonality</b>	Monoclonal
<b>Purity</b>	≥ 90 %
<b>Buffer</b>	0.01M Phosphate Buffered Saline, pH 7.5, containing 0.1% sodium azide
<b>Cross Reactivity</b>	Reactive with 9 aspergillus species including 6 strains of A. fumigatus. A more detailed table of cross-reactivity is available upon request
<b>Specificity</b>	Secreted hyphal antigen from all Aspergillus species tested
<b>Matched Pairs</b>	Self pairs in ELISA
<b>Shelf Life</b>	4 years