

## Product Datasheet

### Anti-Fungal beta-glucan V8201

<b>Product Name</b>	Anti-Fungal beta-glucan
<b>Catalog Number</b>	V8201
<b>Species/Host</b>	Rabbit
<b>Brand</b>	ViroStat
<b>Form/Appearance</b>	Purified from culture medium by protein A chromatography
<b>Concentration</b>	≥ 1 mg/ml
<b>Storage</b>	2-8°C or -20°C. Avoid freeze-thaws. Maximum one freeze-thaw during storage.
<b>Note</b>	Extraction of specimen with 0.4M NaOH required for optimum sensitivity.
<b>Isotype</b>	IgG1
<b>Clonality</b>	Monoclonal
<b>Purity</b>	≥ 90%
<b>Buffer</b>	0.01M Phosphate Buffered Saline, pH 7.5, containing 0.1% sodium azide
<b>Cross Reactivity</b>	Reactive with Candida (4 species), Aspergillus (4 species), Penicillium notatum, Saccharomyces cerevisiae, Pichia pastoris, and Mucor hiemalis. Others not tested.
<b>Specificity</b>	Numerous fungal species
<b>Shelf Life</b>	Four years from date of manufacture