

Product Datasheet

Anti-Fungal beta-glucan V8201

Product Name	Anti-Fungal beta-glucan
Catalog Number	V8201
Description	Purified from culture medium by protein A chromatography
Species/Host	Rabbit
Tested Applications	ELISA
Brand	ViroStat
Form/Appearance	Purified from culture medium by protein A chromatography
Concentration	≥ 1 mg/ml
Storage	+2-8 °C or -20 °C. Avoid freeze-thaws. Maximum one freeze-thaw during storage.
Note	Extraction of specimen with 0.4M NaOH required for optimum sensitivity.
Isotype	IgG1
Clonality	Monoclonal
Purity	≥ 90 %
Buffer	0.01M Phosphate Buffered Saline, pH 7.5, containing 0.1% sodium azide
Cross Reactivity	Reactive with Candida (4 species), Aspergillus (4 species), Penicillium notatum, Saccharomyces cerevisiae, Pichia pastoris, and Mucor hiemalis. Others not tested.
Specificity	Numerous fungal species
Shelf Life	4 years