

Product Datasheet

Anti-Legionella pneumophila V6026

Product Name	Anti-Legionella pneumophila
Catalog Number	V6026
Species/Host	Rabbit
Tested Applications	LF, ELISA, IFA
Brand	ViroStat
Form/Appearance	Purified from culture medium by protein A chromatography
Concentration	≥ 1 mg/ml
Storage	2-8°C or -20°C. Avoid freeze-thaws. Maximum one freeze-thaw during storage.
Isotype	IgG1
Clonality	Monoclonal
Purity	≥ 90%
Buffer	0.01M Phosphate Buffered Saline, pH 7.5, containing 0.1% sodium azide
Cross Reactivity	Non-reactive with human cells
Specificity	Specific to all serogroup 1 strains tested (8)
Shelf Life	Four years from date of manufacture