

## Product Datasheet

### Anti-Legionella pneumophila V6026

<b>Product Name</b>	Anti-Legionella pneumophila
<b>Catalog Number</b>	V6026
<b>Description</b>	Purified from culture medium by protein A chromatography
<b>Species/Host</b>	Rabbit
<b>Tested Applications</b>	ELISA, LF, IFA
<b>Brand</b>	ViroStat
<b>Form/Appearance</b>	Purified from culture medium by protein A chromatography
<b>Concentration</b>	≥ 1 mg/ml
<b>Storage</b>	+2-8 °C or -20 °C. Avoid freeze-thaws. Maximum one freeze-thaw during storage.
<b>Isotype</b>	IgG1
<b>Clonality</b>	Monoclonal
<b>Purity</b>	≥ 90 %
<b>Buffer</b>	0.01M Phosphate Buffered Saline, pH 7.5, containing 0.1% sodium azide
<b>Cross Reactivity</b>	Non-reactive with human cells
<b>Specificity</b>	Specific to all serogroup 1 strains tested (8)
<b>Shelf Life</b>	4 years