

Product Datasheet

Anti-Streptococcus pneumoniae V0494

| | |
|----------------------------|--|
| Product Name | Anti-Streptococcus pneumoniae |
| Catalog Number | V0494 |
| Description | Purified from culture medium by protein A chromatography |
| Species/Host | Rabbit |
| Tested Applications | ELISA, LF, IFA |
| Brand | ViroStat |
| Form/Appearance | Purified from culture medium by protein A chromatography |
| Concentration | ≥ 1 mg/ml |
| Storage | +2-8 °C or -20 °C. Avoid freeze-thaws. Maximum one freeze-thaw during storage. |
| Isotype | IgG1 |
| Clonality | Monoclonal |
| Purity | ≥ 90 % |
| Buffer | 0.01M Phosphate Buffered Saline, pH 7.5, containing 0.1% sodium azide |
| Cross Reactivity | Not determined |
| Specificity | Common cell wall polysaccharide (urinary antigen) of S. pneumoniae |
| Matched Pairs | ELISA pairs (capture-conjugate): V0494-V0495 |
| Shelf Life | 4 years |