

Product Datasheet

Anti-Streptococcus pneumoniae V0495

Product Name	Anti-Streptococcus pneumoniae
Catalog Number	V0495
Description	Purified from culture medium by protein A chromatography
Species/Host	Rabbit
Tested Applications	ELISA, LF, IFA
Brand	ViroStat
Form/Appearance	Purified from culture medium by protein A chromatography
Concentration	≥ 1 mg/ml
Storage	+2-8 °C or -20 °C. Avoid freeze-thaws. Maximum one freeze-thaw during storage.
Isotype	IgG1
Clonality	Monoclonal
Purity	≥ 90 %
Buffer	0.01M Phosphate Buffered Saline, pH 7.5, containing 0.1% sodium azide
Cross Reactivity	Not determined
Specificity	Common cell wall polysaccharide (urinary antigen) of S. pneumoniae
Matched Pairs	ELISA pairs (capture-conjugate): V0494-V0495
Shelf Life	4 years