

# Product Datasheet

## Anti-Streptococcus Group A V0435

Product Name	Anti-Streptococcus Group A
Catalog Number	V0435
Species/Host	Rabbit
Tested Applications	LF, ELISA, IFA
Brand	ViroStat
Form/Appearance	Purified from culture medium by protein A chromatography
Concentration	≥ 1 mg/ml
Storage	2-8°C or -20°C. Avoid freeze-thaws. Maximum one freeze-thaw during storage.
Isotype	IgG1
Clonality	Monoclonal
Purity	≥ 90%
Buffer	0.01M Phosphate Buffered Saline, pH 7.5, containing 0.1% sodium azide
Cross Reactivity	Non-reactive with Strep groups B,C,D,F,G and other bacteria
Specificity	Group A Strep carbohydrate
Matched Pairs	ELISA pairs (capture-conjugate): V0435-V0435, V0435-V1511, V0435-V1501
Shelf Life	Four years from date of manufacture