

# Product Datasheet

## Anti-Streptococcus Group B V1526

Product Name	Anti-Streptococcus Group B
Catalog Number	V1526
Species/Host	Rabbit
Tested Applications	ELISA for detection of extracted group B carbohydrate antigen
Brand	ViroStat
Form/Appearance	Purified from culture medium by protein A chromatography
Concentration	≥ 1 mg/ml
Storage	2-8°C or -20°C. Avoid freeze-thaws. Maximum one freeze-thaw during storage.
Isotype	IgG1
Clonality	Monoclonal
Purity	≥ 90%
Buffer	0.01M Phosphate Buffered Saline, pH 7.5, containing 0.1% sodium azide
Cross Reactivity	Non-reactive with other Strep groups (A, C, D, F & G) and other bacteria
Specificity	Specific to the Group B Strep carbohydrate
Matched Pairs	Self-pairing in ELISA
Shelf Life	Four years from date of manufacture