

Product Datasheet

Anti-Canine Heartworm V3221

Product Name Anti-Canine Heartworm

Catalog Number V3221

Species/Host Goat

Tested Applications ELISA

Brand ViroStat

Form/Appearance Purified IgG fraction prepared by methods developed by Medix

Concentration 5 mg/ml

Storage Short term (<6 months): 2-8°C Long term: -20°C

Isotype IgG

Clonality Polyclonal

Purity ≥95%

Buffer 0.01M Phosphate Buffered Saline, pH 7.2, containing 0.1% sodium azide

Cross Reactivity Antiserum is not absorbed for related filarial Nematodes.
Non-reactive with canine serum

Specificity Soluble & structural antigens of the adult worm

Shelf Life Six years from date of manufacture