

Product Datasheet

Anti-Canine Heartworm FITC V3223

Product Name	Anti-Canine Heartworm FITC
Catalog Number	V3223
Species/Host	Goat
Tested Applications	ELISA, IFA
Brand	ViroStat
Form/Appearance	Purified from antiserum by protein A and covalently coupled with high purity Isomer I of fluorescein isothiocyanate (FITC)
Concentration	5 mg/ml
Storage	Short term (<6 months): 2-8°C Long term: -20°C
Isotype	IgG
Clonality	Polyclonal
Purity	≥95%
Buffer	0.01M Phosphate Buffered Saline, pH 7.2, containing 0.1% sodium azide
Cross Reactivity	Antiserum is not absorbed for related filarial Nematodes. Non-reactive with canine serum
Specificity	Soluble & structural antigens of the adult worm
Shelf Life	Two years from date of manufacture