

Product Datasheet

Anti-Human CD264 Monoclonal Antibody, Azide Free Clone B-T32 853.880.000

Product Name	Anti-Human CD264 Monoclonal Antibody, Azide Free Clone B-T32
Catalog Number	853.880.000
Description	Mouse monoclonal
Species/Host	Mouse
Conjugation	Azide Free
Tested Applications	ELISA
Brand	Diaclone
Form/Appearance	Phosphate-buffered saline. Sterile-filtered through 0.22 µm. Carrier and preservative free
Storage	Stable at +2-8°C for 12 months. For longer storage freeze aliquots.
Isotype	Mouse IgG1 Kappa light chain
Clonality	Monoclonal
Source	Recombinant human TRAIL R4/Fc chimera
Specificity	Recognises TRAIL R4 (DcR2), a 44.2 kDa protein
Hybridoma	Myeloma X63/AG.8653 x Balb/c spleen cells
Immunisation	Recombinant human TRAIL R4/Fc chimera