

# Product Datasheet

## Anti-Papillomavirus V5001

Product Name	Anti-Papillomavirus
Catalog Number	V5001
Description	Purified IgG fraction prepared by methods developed by Medix Biochemica
Species/Host	Goat
Tested Applications	IHC, ELISA, IFA
Brand	ViroStat
Form/Appearance	Purified IgG fraction prepared by methods developed by Medix Biochemica
Concentration	5 mg/ml
Storage	Short term (<6 months): +2-8 °C Long term: -20 °C
Isotype	IgG
Clonality	Polyclonal
Purity	≥ 95 %
Buffer	0.01M Phosphate Buffered Saline, pH 7.2, containing 0.1% sodium azide
Cross Reactivity	Reactive with HPV types 16 & 11 (other types not tested). Reactive in ELISA with Bovine papillomavirus & HPV VLPS. Non-reactive with normal tissue (IHC) & cultured cells (IFA)
Specificity	group
Shelf Life	6 years