

Product Datasheet

Anti-Feline Leukemia virus FITC V7223

| | |
|---------------------|--|
| Product Name | Anti-Feline Leukemia virus FITC |
| Catalog Number | V7223 |
| Species/Host | Goat |
| Tested Applications | IFA, WB immunoassay |
| Brand | ViroStat |
| Form/Appearance | Purified from antiserum by protein A and covalently coupled with high purity Isomer I of fluorescein isothiocyanate (FITC) |
| Concentration | 5 mg/ml |
| Storage | Short term (<6 months): 2-8°C Long term: -20°C |
| Isotype | IgG |
| Clonality | Polyclonal |
| Purity | ≥95% |
| Buffer | 0.01M Phosphate Buffered Saline, pH 7.2, containing 0.1% sodium azide |
| Cross Reactivity | Non-reactive with normal feline serum |
| Specificity | p27 (by western blot) |
| Shelf Life | Two years from date of manufacture |