

## Product Datasheet

### Anti-Herpes simplex virus, types 1&2 FITC V0103

|                            |  |
|----------------------------|--|
| <b>Product Name</b>        | Anti-Herpes simplex virus, types 1&2 FITC  |
| <b>Catalog Number</b>      | V0103  |
| <b>Species/Host</b>        | Rabbit   |
| <b>Tested Applications</b> | IFA, Indirect immunofluorescence   |
| <b>Brand</b>               | ViroStat   |
| <b>Form/Appearance</b>     | Purified from antiserum by protein A and covalently coupled with high purity Isomer I of fluorescein isothiocyanate (FITC) |
| <b>Concentration</b>       | 5 mg/ml  |
| <b>Storage</b>             | Short term (<6 months): 2-8°C Long term: -20°C   |
| <b>Isotype</b>             | IgG  |
| <b>Clonality</b>           | Polyclonal   |
| <b>Purity</b>              | ≥95%   |
| <b>Buffer</b>              | 0.01M Phosphate Buffered Saline, pH 7.2, containing 0.1% sodium azide  |
| <b>Cross Reactivity</b>    | Reactive with types 1 & 2. Non-reactive with HEp-2 cells   |
| <b>Specificity</b>         | ICPs + late structural (virion) antigens   |
| <b>Shelf Life</b>          | Two years from date of manufacture   |