

## Product Datasheet

### Anti-Herpes simplex virus, types 1&2 FITC V0103

<b>Product Name</b>	Anti-Herpes simplex virus, types 1&2 FITC
<b>Catalog Number</b>	V0103
<b>Species/Host</b>	Rabbit
<b>Tested Applications</b>	IFA, Indirect immunofluorescence
<b>Brand</b>	ViroStat
<b>Form/Appearance</b>	Purified from antiserum by protein A and covalently coupled with high purity Isomer I of fluorescein isothiocyanate (FITC)
<b>Concentration</b>	5 mg/ml
<b>Storage</b>	Short term (<6 months): 2-8°C Long term: -20°C
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Purity</b>	≥95%
<b>Buffer</b>	0.01M Phosphate Buffered Saline, pH 7.2, containing 0.1% sodium azide
<b>Cross Reactivity</b>	Reactive with types 1 & 2. Non-reactive with HEp-2 cells
<b>Specificity</b>	ICPs + late structural (virion) antigens
<b>Shelf Life</b>	Two years from date of manufacture