

# Product Datasheet

## Anti-Herpes simplex virus type 2 Biotin V0127

Product Name	Anti-Herpes simplex virus type 2 Biotin
Catalog Number	V0127
Species/Host	Sheep
Tested Applications	ELISA, IFA, WB immunoassay, Indirect immunofluorescence
Brand	ViroStat
Form/Appearance	Purified from antiserum by protein A and covalently coupled with the N-Hydroxysuccinimide ester of biotin under mild conditions to give a high degree of substitution.
Concentration	5 mg/ml
Storage	Short term (<6 months): 2-8°C Long term: -20°C
Isotype	IgG
Clonality	Polyclonal
Purity	≥95%
Buffer	0.01M Phosphate Buffered Saline, pH 7.2, containing 0.1% sodium azide
Cross Reactivity	Non-reactive with HEp-2 cells
Specificity	ICPs + late structural (virion) antigens
Shelf Life	Two years from date of manufacture