

Product Datasheet

Anti-Herpes simplex virus type 2 Biotin V0127

| | |
|----------------------------|--|
| Product Name | Anti-Herpes simplex virus type 2 Biotin |
| Catalog Number | V0127 |
| Description | Purified from antiserum by protein A and covalently coupled with the N-Hydroxysuccinimide ester of biotin under mild conditions to give a high degree of substitution. |
| Species/Host | Sheep |
| Tested Applications | WB, ELISA, IFA, Indirect IFA |
| Brand | ViroStat |
| Form/Appearance | Purified from antiserum by protein A and covalently coupled with the N-Hydroxysuccinimide ester of biotin under mild conditions to give a high degree of substitution. |
| Concentration | 5 mg/ml |
| Storage | Short term (<6 months): +2-8 °C Long term: -20 °C |
| Isotype | IgG |
| Clonality | Polyclonal |
| Purity | ≥ 95 % |
| Buffer | 0.01M Phosphate Buffered Saline, pH 7.2, containing 0.1% sodium azide |
| Cross Reactivity | Non-reactive with HEp-2 cells |
| Specificity | ICPs + late structural (virion) antigens |
| Shelf Life | 2 years |