

# Product Datasheet

## Anti-Herpes simplex virus type 1 FITC V0133

Product Name	Anti-Herpes simplex virus type 1 FITC
Catalog Number	V0133
Species/Host	Goat
Tested Applications	IFA, Indirect immunofluorescence
Brand	ViroStat
Form/Appearance	Purified from antiserum and covalently coupled with high purity Isomer I of fluorescein isothiocyanate (FITC)
Concentration	5 mg/ml
Storage	Short term (<6 months): 2-8°C Long term: -20°C
Isotype	IgG
Clonality	Polyclonal
Purity	≥95%
Buffer	0.01M Phosphate Buffered Saline, pH 7.2, containing 0.1% sodium azide
Cross Reactivity	Reactive with HSV, type 2 by indirect immunofluorescence. Non-reactive with HEp-2 cells
Specificity	ICPs + late structural (virion) antigens
Shelf Life	Two years from date of manufacture