

Product Datasheet

Anti-Herpes simplex virus type 1 Biotin V0137

Product Name	Anti-Herpes simplex virus type 1 Biotin
Catalog Number	V0137
Description	Purified from antiserum and covalently coupled with the N-Hydroxysuccinimide ester of biotin under mild conditions to give a high degree of substitution.
Species/Host	Goat
Tested Applications	ELISA, Indirect IFA
Brand	ViroStat
Form/Appearance	Purified from antiserum and covalently coupled with the N-Hydroxysuccinimide ester of biotin under mild conditions to give a high degree of substitution.
Concentration	5 mg/ml
Storage	Short term (<6 months): +2-8 °C Long term: -20 °C
Isotype	IgG
Clonality	Polyclonal
Purity	≥ 95 %
Buffer	0.01M Phosphate Buffered Saline, pH 7.2, containing 0.1% sodium azide
Cross Reactivity	Reactive with HSV, type 2 by indirect immunofluorescence. Non-reactive with HEp-2 cells
Specificity	ICPs + late structural (virion) antigens
Shelf Life	2 years