

Product Datasheet

Anti-Toxoplasma gondii V0201

Product Name	Anti-Toxoplasma gondii
Catalog Number	V0201
Species/Host	Rabbit
Tested Applications	WB, Indirect immunofluorescence
Brand	ViroStat
Form/Appearance	Purified from antiserum by protein A chromatography
Concentration	5 mg/ml
Storage	Short term (<6 months): 2-8°C Long term: -20°C
Isotype	IgG
Clonality	Polyclonal
Purity	≥95%
Buffer	0.01M Phosphate Buffered Saline, pH 7.2, containing 0.1% sodium azide
Cross Reactivity	Reactive with numerous antigens by Western Blot including SAG1
Specificity	Purified tachyzoites
Shelf Life	Six years from date of manufacture