

# Product Datasheet

## Anti-Epstein-Barr virus (EBNA-1) Biotin V0277

Product Name	Anti-Epstein-Barr virus (EBNA-1) Biotin
Catalog Number	V0277
Species/Host	Goat
Tested Applications	ELISA, ACIF
Brand	ViroStat
Form/Appearance	Purified from antiserum and covalently coupled with the N-Hydroxysuccinimide ester of biotin under mild conditions to give a high degree of substitution.
Concentration	5 mg/ml
Storage	Short term (<6 months): 2-8°C Long term: -20°C
Isotype	IgG
Clonality	Polyclonal
Purity	≥95%
Buffer	0.01M Phosphate Buffered Saline, pH 7.2, containing 0.1% sodium azide
Cross Reactivity	Non-reactive with human lymphoid cells
Specificity	Epstein-Barr nuclear antigen
Shelf Life	Two years from date of manufacture