

Product Datasheet

Anti-Leptospira HRP V0404

Product Name	Anti-Leptospira HRP
Catalog Number	V0404
Species/Host	Rabbit
Tested Applications	ELISA
Brand	ViroStat
Form/Appearance	Purified from antiserum by protein A and covalently coupled to a highly purified preparation of horseradish peroxidase (HRP; RZ>3) by methods developed by Medix Biochemica
Concentration	1 mg/ml
Storage	Short term (<6 months): 2-8°C Long term: -20°C
Isotype	IgG
Clonality	Polyclonal
Purity	≥95%
Buffer	0.01M Phosphate Buffered Saline, pH 7.2, containing bovine serum albumin (10 mg/ml) as a stabilizer and preserved with 0.002% thimerosal
Cross Reactivity	Reactive in ELISA with L. biflexa, L. australis, L. canicola, L. grippotyphosa, L. icterohaemorrhagiae, L. Pomona, L. serjroe. Antiserum has not been absorbed and may react with related microorganisms
Specificity	All antigens
Shelf Life	Two years from date of manufacture