

Product Datasheet

Anti-Leptospira HRP V0404

Product Name	Anti-Leptospira HRP
Catalog Number	V0404
Description	Purified from antiserum by protein A and covalently coupled to a highly purified preparation of horseradish peroxidase (HRP; RZ>3) by methods developed by Medix Biochemica
Species/Host	Rabbit
Tested Applications	ELISA
Brand	ViroStat
Form/Appearance	Purified from antiserum by protein A and covalently coupled to a highly purified preparation of horseradish peroxidase (HRP; RZ>3) by methods developed by Medix Biochemica
Concentration	1 mg/ml
Storage	Short term (<6 months): +2-8 °C Long term: -20 °C
Isotype	IgG
Clonality	Polyclonal
Purity	≥ 95 %
Buffer	0.01M Phosphate Buffered Saline, pH 7.2, containing bovine serum albumin (10 mg/ml) as a stabilizer and preserved with 0.002% thimerosal
Cross Reactivity	Reactive in ELISA with <i>L. biflexa</i> , <i>L. australis</i> , <i>L. canicola</i> , <i>L. grippotyphosa</i> , <i>L. icterohaemorrhagiae</i> , <i>L. Pomona</i> , <i>L. serjroe</i> . Antiserum has not been absorbed and may react with related microorganisms
Specificity	All antigens
Shelf Life	2 years