

Product Datasheet

Anti-Streptococcus pneumoniae FITC V0453

Product Name	Anti-Streptococcus pneumoniae FITC
Catalog Number	V0453
Description	Purified from antiserum by protein A and covalently coupled with high purity Isomer I of fluorescein isothiocyanate (FITC)
Species/Host	Rabbit
Tested Applications	IFA
Brand	ViroStat
Form/Appearance	Purified from antiserum by protein A and covalently coupled with high purity Isomer I of fluorescein isothiocyanate (FITC)
Concentration	5 mg/ml
Storage	Short term (<6 months): +2-8 °C Long term: -20 °C
Isotype	IgG
Clonality	Polyclonal
Purity	≥ 95 %
Buffer	0.01M Phosphate Buffered Saline, pH 7.2, containing 0.1% sodium azide
Cross Reactivity	Antiserum is not absorbed and may react with related microorganisms. Reactive with types 3, 4, 6, 7, 9, 14, 18, 19 & 23
Specificity	Types 3, 4, 6, 7, 9, 14, 18, 19, & 23
Shelf Life	2 years