

Product Datasheet

Anti-Streptococcus pneumoniae HRP V0454

| | |
|---------------------|---|
| Product Name | Anti-Streptococcus pneumoniae HRP |
| Catalog Number | V0454 |
| Species/Host | Rabbit |
| Tested Applications | ELISA |
| Brand | ViroStat |
| Form/Appearance | Purified from antiserum by protein A and covalently coupled to a highly purified preparation of horseradish peroxidase (HRP; RZ>3) by methods developed by Medix Biochemica |
| Concentration | 1 mg/ml |
| Storage | Short term (<6 months): 2-8°C Long term: -20°C |
| Isotype | IgG |
| Clonality | Polyclonal |
| Purity | ≥95% |
| Buffer | 0.01M Phosphate Buffered Saline, pH 7.2, containing bovine serum albumin (10 mg/ml) as a stabilizer and preserved with 0.002% thimerosal |
| Cross Reactivity | Antiserum is not absorbed and may react with related microorganisms. Reactive with types 3, 4, 6, 7, 9, 14, 18, 19 & 23 |
| Specificity | Types 3, 4, 6, 7, 9, 14, 18, 19, & 23 |
| Shelf Life | Two years from date of manufacture |