

Product Datasheet

Anti-Streptococcus pneumoniae HRP V0454

Product Name	Anti-Streptococcus pneumoniae HRP
Catalog Number	V0454
Description	Purified from antiserum by protein A and covalently coupled to a highly purified preparation of horseradish peroxidase (HRP; RZ>3) by methods developed by Medix Biochemica
Species/Host	Rabbit
Tested Applications	ELISA
Brand	ViroStat
Form/Appearance	Purified from antiserum by protein A and covalently coupled to a highly purified preparation of horseradish peroxidase (HRP; RZ>3) by methods developed by Medix Biochemica
Concentration	1 mg/ml
Storage	Short term (<6 months): +2-8 °C Long term: -20 °C
Isotype	IgG
Clonality	Polyclonal
Purity	≥ 95 %
Buffer	0.01M Phosphate Buffered Saline, pH 7.2, containing bovine serum albumin (10 mg/ml) as a stabilizer and preserved with 0.002% thimerosal
Cross Reactivity	Antiserum is not absorbed and may react with related microorganisms. Reactive with types 3, 4, 6, 7, 9, 14, 18, 19 & 23
Specificity	Types 3, 4, 6, 7, 9, 14, 18, 19, & 23
Shelf Life	2 years