

Product Datasheet

Anti-Rotavirus V0501

Product Name	Anti-Rotavirus
Catalog Number	V0501
Species/Host	Goat
Tested Applications	Indirect immunofluorescence
Brand	ViroStat
Form/Appearance	Purified IgG fraction prepared by methods developed by Medix Biochemica
Concentration	5 mg/ml
Storage	Short term (<6 months): 2-8°C Long term: -20°C
Isotype	IgG
Clonality	Polyclonal
Purity	≥95%
Buffer	0.01M Phosphate Buffered Saline, pH 7.2, containing 0.1% sodium azide
Cross Reactivity	>90% with human rotaviruses. Non-reactive with HEp-2 cells and WI-38 cells by indirect immunofluorescence
Specificity	ICPs + late structural (virion) antigens
Shelf Life	Six years from date of manufacture