

# Product Datasheet

## Anti-Rotavirus FITC V0503

Product Name	Anti-Rotavirus FITC
Catalog Number	V0503
Description	Purified from antiserum and covalently coupled with high purity Isomer I of fluorescein isothiocyanate (FITC)
Species/Host	Goat
Tested Applications	IFA, Indirect IFA
Brand	ViroStat
Form/Appearance	Purified from antiserum and covalently coupled with high purity Isomer I of fluorescein isothiocyanate (FITC)
Concentration	5 mg/ml
Storage	Short term (<6 months): +2-8 °C Long term: -20 °C
Isotype	IgG
Clonality	Polyclonal
Purity	≥ 95 %
Buffer	0.01M Phosphate Buffered Saline, pH 7.2, containing 0.1% sodium azide
Cross Reactivity	>90% with human rotaviruses. Non-reactive with HEp-2 cells and WI-38 cells by indirect immunofluorescence
Specificity	ICPs + late structural (virion) antigens
Shelf Life	2 years