

Product Datasheet

Anti-Respiratory Syncytial virus HRP
V0604

Product Name	Anti-Respiratory Syncytial virus HRP
Catalog Number	V0604
Species/Host	Goat
Tested Applications	ELISA
Brand	ViroStat
Form/Appearance	Purified from antiserum and covalently coupled to a highly purified preparation of horseradish peroxidase (HRP; RZ>3) by methods developed by Medix Biochemica
Concentration	1 mg/ml
Storage	Short term (<6 months): 2-8°C Long term: -20°C
Isotype	IgG
Clonality	Polyclonal
Purity	≥95%
Buffer	0.01M Phosphate Buffered Saline, pH 7.2, containing bovine serum albumin (10 mg/ml) as a stabilizer and preserved with 0.002% thimerosal
Cross Reactivity	Non-reactive with Para 1-3, Influenza A & B, & Adenovirus. Non-reactive with HEp-2 and WI-38 cells
Specificity	All viral antigens
Shelf Life	Two years from date of manufacture