

Product Datasheet

Anti-Respiratory Syncytial virus Biotin V0607

Product Name	Anti-Respiratory Syncytial virus Biotin
Catalog Number	V0607
Species/Host	Goat
Tested Applications	ELISA, IFA, immunohistochemistry
Brand	ViroStat
Form/Appearance	Purified from antiserum and covalently coupled with the N-Hydroxysuccinimide ester of biotin under mild conditions to give a high degree of substitution.
Concentration	5 mg/ml
Storage	Short term (<6 months): 2-8°C Long term: -20°C
Isotype	IgG
Clonality	Polyclonal
Purity	≥95%
Buffer	0.01M Phosphate Buffered Saline, pH 7.2, containing 0.1% sodium azide
Cross Reactivity	Non-reactive with Para 1-3, Influenza A & B, & Adenovirus by IFA. Non-reactive with HEp-2 and WI-38 cells
Specificity	All viral antigens
Shelf Life	Two years from date of manufacture