

Product Datasheet

Anti-Salmonella V0701

Product Name	Anti-Salmonella
Catalog Number	V0701
Description	Purified from antiserum by protein A chromatography
Species/Host	Rabbit
Tested Applications	Indirect IFA
Brand	ViroStat
Form/Appearance	Purified from antiserum by protein A chromatography
Concentration	5 mg/ml
Storage	Short term (<6 months): +2-8 °C Long term: -20 °C
Isotype	IgG
Clonality	Polyclonal
Purity	≥ 95 %
Buffer	0.01M Phosphate Buffered Saline, pH 7.2, containing 0.1% sodium azide
Cross Reactivity	Antiserum is unabsorbed and does cross react with related Enterobacteriaceae
Specificity	Polyvalent for Salmonella "O" & "H" antigens
Shelf Life	6 years