

# Product Datasheet

## Anti-Shigella V0901

|                     |  |
|---------------------|--|
| Product Name        | Anti-Shigella  |
| Catalog Number      | V0901  |
| Species/Host        | Rabbit   |
| Tested Applications | ELISA, IFA   |
| Brand               | ViroStat   |
| Form/Appearance     | Purified from antiserum by protein A chromatography                      |
| Concentration       | 5 mg/ml  |
| Storage             | Short term (<6 months): 2-8°C Long term: -20°C                           |
| Isotype             | IgG  |
| Clonality           | Polyclonal   |
| Purity              | ≥95%   |
| Buffer              | 0.01M Phosphate Buffered Saline, pH 7.2, containing 0.1% sodium azide    |
| Cross Reactivity    | Antiserum is not absorbed and does react with related Enterobacteriaceae |
| Specificity         | Polyvalent for Shigella antigens   |
| Shelf Life          | Six years from date of manufacture                                       |