

# Product Datasheet

## Anti-Shigella V0901

Product Name	Anti-Shigella
Catalog Number	V0901
Species/Host	Rabbit
Tested Applications	ELISA, IFA
Brand	ViroStat
Form/Appearance	Purified from antiserum by protein A chromatography
Concentration	5 mg/ml
Storage	Short term (<6 months): 2-8°C Long term: -20°C
Isotype	IgG
Clonality	Polyclonal
Purity	≥95%
Buffer	0.01M Phosphate Buffered Saline, pH 7.2, containing 0.1% sodium azide
Cross Reactivity	Antiserum is not absorbed and does react with related Enterobacteriaceae
Specificity	Polyvalent for Shigella antigens
Shelf Life	Six years from date of manufacture