

Product Datasheet

Anti-E. coli V1001

Product Name	Anti-E. coli
Catalog Number	V1001
Species/Host	Rabbit
Tested Applications	ELISA
Brand	ViroStat
Form/Appearance	Purified from antiserum by protein A chromatography
Concentration	5 mg/ml
Storage	Short term (<6 months): 2-8°C Long term: -20°C
Isotype	IgG
Clonality	Polyclonal
Purity	≥95%
Buffer	0.01M Phosphate Buffered Saline, pH 7.2, containing 0.1% sodium azide
Cross Reactivity	Antiserum is not absorbed and does react with related Enterobacteriaceae
Specificity	All “O” and “K” antigenic serotypes of Escherichia coli
Shelf Life	Six years from date of manufacture