

# Product Datasheet

## Anti-E. coli FITC V1003

Product Name	Anti-E. coli FITC
Catalog Number	V1003
Species/Host	Rabbit
Tested Applications	ELISA, IFA
Brand	ViroStat
Form/Appearance	Purified from antiserum by protein A and covalently coupled with high purity Isomer I of fluorescein isothiocyanate (FITC)
Concentration	5 mg/ml
Storage	Short term (<6 months): 2-8°C Long term: -20°C
Isotype	IgG
Clonality	Polyclonal
Purity	≥95%
Buffer	0.01M Phosphate Buffered Saline, pH 7.2, containing 0.1% sodium azide
Cross Reactivity	Antiserum is not absorbed and does react with related Enterobacteriaceae
Specificity	All “O” and “K” antigenic serotypes of Escherichia coli
Shelf Life	Two years from date of manufacture