

Product Datasheet

Anti-E. coli V1091

Product Name	Anti-E. coli
Catalog Number	V1091
Species/Host	Goat
Tested Applications	Indirect IFA
Brand	ViroStat
Form/Appearance	This product consists of the purified IgG fraction prepared by methods developed by Medix Biochemica
Concentration	5 mg/ml
Storage	Short term (<6 months): 2-8°C Long term: -20°C
Isotype	IgG
Clonality	Polyclonal
Purity	≥95%
Buffer	0.01M Phosphate Buffered Saline, pH 7.2, containing 0.1% sodium azide
Cross Reactivity	Antiserum is unabsorbed and does cross react with related Enterobacteriaceae
Specificity	Many “O” and “K” antigenic serotypes of Escherichia coli
Shelf Life	Six years from date of manufacture