

Product Datasheet

Anti-E. coli V1091

Product Name	Anti-E. coli
Catalog Number	V1091
Description	This product consists of the purified IgG fraction prepared by methods developed by Medix Biochemica
Species/Host	Goat
Tested Applications	Indirect IFA
Brand	ViroStat
Form/Appearance	This product consists of the purified IgG fraction prepared by methods developed by Medix Biochemica
Concentration	5 mg/ml
Storage	Short term (<6 months): +2-8 °C Long term: -20 °C
Isotype	IgG
Clonality	Polyclonal
Purity	≥ 95 %
Buffer	0.01M Phosphate Buffered Saline, pH 7.2, containing 0.1% sodium azide
Cross Reactivity	Antiserum is unabsorbed and does cross react with related Enterobacteriaceae
Specificity	Many "O" and "K" antigenic serotypes of Escherichia coli
Shelf Life	6 years