

Product Datasheet

Anti-E. coli FITC
V1093

Product Name	Anti-E. coli FITC
Catalog Number	V1093
Species/Host	Goat
Tested Applications	Indirect IFA
Brand	ViroStat
Form/Appearance	Purified from antiserum by protein A and covalently coupled with high purity Isomer I of fluorescein isothiocyanate (FITC)
Concentration	5 mg/ml
Storage	Short term (<6 months): 2-8°C Long term: -20°C
Isotype	IgG
Clonality	Polyclonal
Purity	≥95%
Buffer	0.01M Phosphate Buffered Saline, pH 7.2, containing 0.1% sodium azide
Cross Reactivity	Antiserum is unabsorbed and does cross react with related Enterobacteriaceae
Specificity	Many “O” and “K” antigenic serotypes of Escherichia coli
Shelf Life	Two years from date of manufacture