

# Product Datasheet

## Anti-E. coli HRP V1094

Product Name	Anti-E. coli HRP
Catalog Number	V1094
Species/Host	Goat
Tested Applications	ELISA
Brand	ViroStat
Form/Appearance	Purified from antiserum and covalently coupled to a highly purified preparation of horseradish peroxidase (HRP; RZ>3) by methods developed by Medix Biochemica
Concentration	1 mg/ml
Storage	Short term (<6 months): 2-8°C Long term: -20°C
Isotype	IgG
Clonality	Polyclonal
Purity	≥95%
Buffer	0.01M Phosphate Buffered Saline, pH 7.2, containing bovine serum albumin (10 mg/ml) as a stabilizer and preserved with 0.002% thimerosal
Cross Reactivity	Antiserum is unabsorbed and does cross react with related Enterobacteriaceae
Specificity	Many “O” and “K” antigenic serotypes of Escherichia coli
Shelf Life	Two years from date of manufacture