

Product Datasheet

Anti-Parainfluenza virus type 3 FITC V1203

Product Name	Anti-Parainfluenza virus type 3 FITC
Catalog Number	V1203
Description	Purified from antiserum by protein A and covalently coupled with high purity Isomer I of fluorescein isothiocyanate (FITC)
Species/Host	Goat
Tested Applications	IFA, Indirect IFA
Brand	ViroStat
Form/Appearance	Purified from antiserum by protein A and covalently coupled with high purity Isomer I of fluorescein isothiocyanate (FITC)
Concentration	5 mg/ml
Storage	Short term (<6 months): +2-8 °C Long term: -20 °C
Isotype	IgG
Clonality	Polyclonal
Purity	≥ 95 %
Buffer	0.01M Phosphate Buffered Saline, pH 7.2, containing 0.1% sodium azide
Cross Reactivity	Reactive with type 2. Minimal reactivity with type 1. Some cross-reactivity with bovine para-3 (= 1:400) & canine para [1,000]. Non-reactive with HEp-2 cells by indirect immunofluorescence
Specificity	All structural antigens
Shelf Life	2 years