

# Product Datasheet

## Anti-Parainfluenza virus type 1 V1211

Product Name	Anti-Parainfluenza virus type 1
Catalog Number	V1211
Species/Host	Goat
Tested Applications	Indirect immunofluorescence, Indirect IFA
Brand	ViroStat
Form/Appearance	Purified IgG fraction prepared by methods developed by Medix Biochemica
Concentration	5 mg/ml
Storage	Short term (<6 months): 2-8°C Long term: -20°C
Isotype	IgG
Clonality	Polyclonal
Purity	≥95%
Buffer	0.01M Phosphate Buffered Saline, pH 7.2, containing 0.1% sodium azide
Cross Reactivity	Non-reactive with Para types 2 & 3 by indirect immunofluorescence. Non-reactive with Hep-2 cells or monkey kidney cells by indirect IFA
Specificity	All structural antigens
Shelf Life	Two years from date of manufacture