

Product Datasheet

Anti-Human IL-23 Monoclonal Antibody, Azide Free Clone
B-F43
853.930

Product Name	Anti-Human IL-23 Monoclonal Antibody, Azide Free Clone B-F43
Catalog Number	853.930
Description	Mouse monoclonal antibody. Source: Recombinant antibody from the original clone
Species/Host	Mouse
Conjugation	Azide Free
Tested Applications	ELISA
Immunogen	Recombinant human IL-23 (sub-unit p19)
Brand	Diaclone
Form/Appearance	Liquid
Concentration	1.0 mg/ml (+/- 10 %)
Storage	Stable at +2-8 °C for 12 months. For longer storage freeze aliquots.
Isotype	Mouse IgG1 Kappa light chain
Clonality	Monoclonal
Source	Recombinant antibody from the original clone
Buffer	Phosphate-buffered saline. Sterile-filtered through 0.22µm. Carrier and preservative free.
Specificity	Recognizes both recombinant and natural human IL-23, p19 sub-unit