

Product Datasheet

Anti-Parainfluenza virus type 1 HRP V1214

Product Name	Anti-Parainfluenza virus type 1 HRP
Catalog Number	V1214
Description	Purified from antiserum and covalently coupled to a highly purified preparation of horseradish peroxidase (HRP; RZ>3) by methods developed by Medix Biochemica
Species/Host	Goat
Tested Applications	Indirect IFA
Brand	ViroStat
Form/Appearance	Purified from antiserum and covalently coupled to a highly purified preparation of horseradish peroxidase (HRP; RZ>3) by methods developed by Medix Biochemica
Concentration	1 mg/ml
Storage	Short term (<6 months): +2-8 °C Long term: -20 °C
Isotype	IgG
Clonality	Polyclonal
Purity	≥ 95 %
Buffer	0.01M Phosphate Buffered Saline, pH 7.2, containing bovine serum albumin (10 mg/ml) as a stabilizer and preserved with 0.002% thimerosal
Cross Reactivity	Non-reactive with Para types 2 & 3. Non-reactive with HEp-2 cells or monkey kidney cells
Specificity	All structural antigens
Shelf Life	2 years