

Product Datasheet

Anti-Parainfluenza virus type 1 Biotin V1217

Product Name	Anti-Parainfluenza virus type 1 Biotin
Catalog Number	V1217
Description	Purified from antiserum and covalently coupled with the N-Hydroxysuccinimide ester of biotin under mild conditions to give a high degree of substitution.
Species/Host	Goat
Tested Applications	ELISA, Indirect IFA
Brand	ViroStat
Form/Appearance	Purified from antiserum and covalently coupled with the N-Hydroxysuccinimide ester of biotin under mild conditions to give a high degree of substitution.
Concentration	5 mg/ml
Storage	Short term (<6 months): +2-8 °C Long term: -20 °C
Isotype	IgG
Clonality	Polyclonal
Purity	≥ 95 %
Buffer	0.01M Phosphate Buffered Saline, pH 7.2, containing 0.1% sodium azide
Cross Reactivity	Non-reactive with Para types 2 & 3 by indirect immunofluorescence. Non-reactive with Hep-2 cells or monkey kidney cells by indirect IFA
Specificity	All structural antigens
Shelf Life	2 years